

JAW CRUSHER



Brand	: ESSA
Model	: JC1000
Jaw Inlet	: 130 mm x 100 mm
Max. Feed Size	: 85 mm
Usual Product Size	: 6 mm
Max. Product Size	: 25 mm
Throughput at Usual Product Size	: 250 kg/h
Motor Power	: 2.2 kW, 3 phase
Power Supply	: AC 380V, 3 phase + earth
Toggle Speed	: 350 stokes per minute
Dimension	: 1115 x 485 x 1180 mm
Weight	: 250 kg