

GRINDING MILL



Vessel includes vessel body, impact ring, and impact block. The sample places between them, the sample diameter is most greatly 5mm



Brand	:	ZHONGHE
Model	:	ZHM-1-T
Input Sample	:	≤ 10 mm
Output Sample	:	≤ 80 μm
Sample Capacity	:	30 - 80 grams
Grinding Time	:	0 - 999 s
Power Supply	:	AC 380 V, 50 Hz, 3 φ
Rated Load	:	1.5 Kw
Rotate Speed	:	960/1400 r/min
Vessel Amount	:	1
Vessel Cubage	:	100 cc
Dimension	:	650 x 600 x 1000 mm
Weight	:	220 Kg
Ambient Temperature:	:	10 - 40 °C
Vessel Material	:	Manganese Steel